

Engineering Conferences International

ECI Digital Archives

Single-Use Technologies V: Building The Future

Proceedings

3-21-2022

Measuring the liquid volume in a flexible bioprocess container with Single Use Pressure Sensors

Nick Troise

James Furey

Follow this and additional works at: https://dc.engconfintl.org/sut_v

Measuring the Liquid Volume in a Flexible Bioprocess Container with Single Use Pressure Sensors



Nick Troise
R&D Engineer II
March, 2022

Phone: +1.609.799.2299
Email: nick.troise@PendoTECH.com

© Copyright 2022



Agenda

- Introduction
- Highlight Technology
- Application Overview
- Background on Pascal's Principle and Hydrostatic Pressure
- Review Validation Experiment
- Analyze Results
- Discuss Key Takeaways
- Areas of Improvement
- Conclusion and Next Steps

Nick Troise


R&D Engineer II

- Graduated from Rowan University in Biomedical Engineering with a dual minor in Math and Computer Science
- Started working with PendoTECH in 2017 as an intern and started full time in 2019
- Responsibilities are focused on new product design/development, and product validation
- Represent PendoTECH on the Bioprocess System's Alliance (BPSA) X-ray Irradiation Committee
- Fun Facts:
 - Avid fan of both Chess and MMA
 - Lost about 70 lbs (32 kg) in last 10 months



Single Use Port Plate Pressure Sensor

- Designed for integration with flexible bioprocess containers (BPC)
 - Buffer/media storage bags, mixing bags, bioreactors, etc.
 - Interfaces with port welded onto BPC
- Single Use Pressure Sensor
 - Liquid or gas applications
 - No calibration required (plug and play)
 - Low cost and can be disposed with rest of assembly
- Key Features
 - Same MEMS-HAP sensing technology as inline pressure sensors
 - Qualified for use up to 15 psig/1.03bar (equivalent to a 34ft/10m tall bag)
 - Fluid contact materials all USP Class VI and ADCF
 - Validated by customers to weld to various bag materials
 - Gamma Irradiation Compatible



Key Applications

- Safety-
 - Preventing over-pressurization of a flexible BPC
- Performance
 - Monitoring pressure in a SU Bioreactor
- Leak Testing
 - Pre-use integrity test via pressure decay
- **Volume Measurement by hydrostatic pressure**



Application Overview

- Measuring volume of a flexible BPC
 - Important for monitoring levels of a buffer or media storage
 - Alert operators when bags are running low or need to be switched
- Current Solutions:
 - Load Cells
 - Very Accurate and high resolution
 - Expensive, can be difficult to integrate, and requires calibration
 - Level Sensors
 - Several different styles: Conductive, ultrasonic, capacitance, etc.
 - Affordable and accurate, but can be affected by factors in the e
 - Generally pretty invasive



Using Pascal's Principle to measure liquid volume

- Define mathematical equation for hydrostatic pressure
- *In a static environment, the depth of a liquid generates a force that is directly proportional to the height of the liquid*
 - Independent of the volume, shape, total mass, or area of fluid
- Given a known pressure and fluid density, can solve for Δh
- With Δh and cross section area (A) can determine the volume of liquid


$$\Delta P = \rho g(\Delta h)$$

ΔP = hydrostatic pressure

ρ = fluid density

g = acceleration due to gravity,

Δh = height of the liquid.





Important Requirements

- No external or head pressure in the system
 - Or at least must be quantified so it can be offset before calculating Δh
- Locational Accuracy of the Sensor
 - Precise location must be known, ideally directly on bottom of bag
 - Otherwise the pressure will not correlate to the true liquid level
- High Resolution Pressure Reading
 - Water column of 1 in (in Aq) = 0.036 psi (2.5 mbar)
 - Must be able to detect very small changes in pressure for accurate level measurements



Validation Experiment

- PendoTECH Port Plate installed on the bottom of a 200L 2D BPC
- Supported by custom 200L tote with cut-out for port plate sensor
- Monitor: PMAT2HR
 - Specializes in low pressure applications
 - Zero offset stability = $\pm 0.003\text{psi}$ (0.21mbar)




Validation Experiment (cont.)

- Filled container dyed water to make physical measurements of liquid level with caliper
 - Slits cut into side to easily identify meniscus of liquid level
 - Measured every 0.20psi or 14mbar (4-6 in Aq) until BPC was full
- Calculated volume using equation for volume of a cylinder: $V=\pi r^2 h$
 - No direct volume reference
 - Graduation marks on tote used as a sanity/reasonableness check



Results Analysis


Bag Pressure Measured (psi)	Theoretical Liquid Level (in Aq)*	Measured Liquid Level (in Aq)*	Percent Error	Calculated Volume (L)
0.305	8.45	8.54	1.07%	53.2
0.451	12.50	12.56	0.48%	78.2
0.620	17.18	17.19	0.06%	107.1
0.769	21.31	21.27	0.19%	132.5
0.920	25.49	25.88	1.53%	161.2
1.13	31.31	31.54	0.73%	196.5



- Theoretical and measured liquid levels extremely similar
 - Differed by an average <1%
- Greatest error was 1.53% at 0.920 psi
 - Within sensor's accuracy spec of $\pm 2\%$ of reading
- Relationship between liquid level and bag pressure was perfectly linear as expected
- Qualitative comparison of volume using graduation marks on tote
 - Bag filled to 160L mark, calculated volume was 166L, $\approx 4\%$ error

Key Observations

- Hydrostatic pressure can be used to measure liquid level very accurately under the right conditions
 - Sensor location precisely determined
 - High resolution pressure reading
- Correlation to total liquid volume is highly dependent on bag geometry
 - Calculation typically assumes an ideal or perfect shape
 - Folds/creases caused by deployment can significantly impact volume calculation
 - Cross Sectional Area key contributor to volume accuracy
 - Greater ratio of area to height, increases potential inaccuracy
 - Even with poor deployment/bag geometry, can still be used for relative volume measurements
 - E.g. Gas Gauge
 - Under ideal settings, can measure volume very accurately
 - Comparable performance to load cell and other technologies



Potential Improvements

- Use a 3D BPC with multiple sensor ports
 - Larger sample size = more data for repeatability/consistency
 - Easier to work with/deploy for volume calculations
- Reliable reference for volume measurement
 - Measuring liquid level is good step, but volume accuracy is the ultimate goal
 - Requires an accurate control to evaluate
 - Calibrated scale below BPC, or flow meter on fill line to measure actual volume



Conclusion and Next Steps

- Following Pascal's Principle, a single use port plate pressure sensor can accurately measure the liquid level in a flexible bioprocess container
 - With correct setup this can be used to accurately measure the volume in a BPC in a cost-effective, simple, and minimally invasive manner
- Complete a re-vamped application validation experiment
 - Implement improvements previously discussed to quantitatively assess volume accuracy
 - Can perform with a gamma irradiated sensor to further challenge gamma compatibility in setting that reflects actual customer usage
- Partner with end-users for additional data and feedback



Thank you