

Fall 11-2-2015

Inline spiking for viral clearance validation of continuous processes

Herb Lutz
EMD Millipore

Follow this and additional works at: http://dc.engconfintl.org/biomanufact_ii



Part of the [Biomedical Engineering and Bioengineering Commons](#)

Recommended Citation

Herb Lutz, "Inline spiking for viral clearance validation of continuous processes" in "Integrated Continuous Biomanufacturing II", Chetan Goudar, Amgen Inc. Suzanne Farid, University College London Christopher Hwang, Genzyme-Sanofi Karol Lacki, Novo Nordisk Eds, ECI Symposium Series, (2015). http://dc.engconfintl.org/biomanufact_ii/97

This Conference Proceeding is brought to you for free and open access by the Proceedings at ECI Digital Archives. It has been accepted for inclusion in Integrated Continuous Biomanufacturing II by an authorized administrator of ECI Digital Archives. For more information, please contact franco@bepress.com.

Inline Spiking Method for Validating Continuous Virus Filtration

Herb Lutz (principal consulting engineer BSN process solutions)

Introduction

Conventional virus filtration (VF) processes involve batch processing a uniform feed at uniform conditions (e.g. constant pressure). Continuous VF processes can involve an elution protein peak with a salt gradient under non-uniform pressure where standard spiking methodologies do not represent the manufacturing process as required by the ICH Q5A virus validation guidelines.

This paper describes the development, qualification, and implementation of an inline spiking methodology that is suitable for validating continuous processes.

Standard spiking

Batch prefiltration

Flush prefilter with water and buffer. Perform prefiltration at constant pressure. Filter at typical process times of ~2hrs to ensure comparable residence times for plugging agent removal



Batch virus filtration

Add virus spike to batch prefiltered product. Filter spiked solution with 0.2 or 0.45mm membranes. Perform virus filtration, sample filtrate weight at regular time points

Inline spiking

Batch prefilter-> spike-> virus filter

Flush filters with water and buffer. Add 0.2 mm filtered protein solution to the feed vessel and 0.2 mm filtered spiking solution to the feed syringe pump. Pressurize feed vessel, open valves. Start syringe pump to continuously inject spike, mix with refiltered feed, sample, and feed to virus filter. Measure filtrate weight vs time & calculate flow rate. As the virus filter plugs and flow drops, slow down the syringe pump to maintain % spiking ratio




Inline Qualification with salt spike

Salt spike improves sensitivity of concentration measurements: 2x difference in salt concentration (conductivity) corresponds to a 0.3 difference in log₁₀ scales.

Run #	% Spike flow to filtrate flow	Spike (M)	Predicted Filtrate/Sample Syringe (M)	Sample Syringe (M)	Filtrate (M)
1	5%	0.68	0.028	0.024	0.028
2	5%	0.68	0.029	0.025	0.029
3	5%	0.68	0.026	0.024	0.026
4	10%	0.68	0.070	0.063	0.064
5	10%	0.68	0.069	0.057	0.064
6	10%	0.68	0.068	0.064	0.065

Run to run variability was <6%. Sample syringe titer can be used as feed pool titer. Sample syringe and filtrate within 14% (0.07 log₁₀ units or ~0 LRV). The predicted filtrate/sample syringe titer within 9% (Hold control ~0 LRV). Well mixed solution so sample syringe is representative of cross section. Mass balance closes & holdup volumes are negligible

Inline qualification with prefiltration mAb2 and Viresolve™ Pro spiked with 5% MMV or 5% xMuLV



MMV plugging due to protein. xMuLV plugging due to protein+spike. 2-4x capacity increases with inline spiking (even with xMuLV spike)


Run Description	Spike Titer log(TCID ₅₀ /mL)	Predicted Sample Titer Log(TCID ₅₀ /mL)	Actual Sample Titer Log(TCID ₅₀ /mL)	LRV
MMV-inline	6.68	5.38	5.30	≥4.42
MMV-standard	--	--	5.68	≥4.80
xMuLV-inline	6.10	4.80	4.60	≥3.42
xMuLV-standard	--	--	4.48	≥3.30

Inline spiking shows comparable LRV to standard spiking. Inline spiking shows improved filtrate throughput


Linked processing

Sequential unit operations are started before the previous steps are completed.

Virus filter or other clearance step sees variable operating feed parameters




ICH Q5A (1996) Measure virus removal using scaled down models that mirror large scale production to ensure that the virus removal results can be translated to the biotherapeutic product.



Linked processing data

Inline spiking enables linked process validation study


>5 LRV consistent with batch processes Feed & sample loads match (no inactivation) Repeatable



Run 1	Titer (LogTCID ₅₀ /mL)	Vol (mL)	Viral Load (Log TCID)	
ISS (feed syr.)	6.88	5.1	7.86	Avg 7.85
IHC (sample syr.)	5.30	112	7.83	
Pool	-0.62	112	1.91	

Run 2	Titer (LogTCID ₅₀ /mL)	Vol (mL)	Viral Load (Log TCID)	
ISS (feed syr.)	6.55	3.1	7.52	Avg 7.47
IHC (sample syr.)	4.93	104	7.42	
Pool	≤ -0.62	104	≤ 1.88	

Data courtesy Megan McClure, Viveka Raol, Amgen



MuLV elution from upstream CEX a) pH5, b) pH6

Summary

Inline spiking methodology:

- Improves filter throughputs
- More representative of manufacturing inline operation
- Demonstrates virus LRV comparable to standard spiking
- Consistent
- Relatively easy to operate
- Useful where proteins form plugging foulants
- Facilitates validation of linked processes

References

- CBER/FDA, "Points to consider in the manufacture and testing of monoclonal antibody products for human use", 1997
- CBER/FDA, "Draft of points to consider in the characterization of cell lines used to produce biologicals", 1993
- Connell-Crowley, L., et al., "Cation exchange chromatography provides effective retrovirus clearance for antibody purification processes", *Biotechnol. Bioeng.*, 109(2012)157-165
- CPMP/EMA, "Note for guidance on virus validation studies: the design, contribution, and interpretation of studies validating the inactivation and removal of viruses", 1996
- ICH, "Viral safety evaluation of biotechnology products", 1997
- Lutz, et al., *Biotechnol. Progress*, 27 (2011) 121-128
- PDA, Technical Report No. 47 Preparation of Virus Spikes Used for Virus Clearance Studies, ISBN: 978-0-939459-27-8, 2010, pp. 199-208.

Acknowledgments

Eva Gefroh, Megan McClure (Just Pharm), Brian Hubbard, Viveka Raol, Joseph Shultz (Amgen), Gisela Ramsey, Joe Parrella, Paul Genest (EMD Millipore), Wayne Chang (Boehringer Ingelheim), Tamas Blandl (MSD), John Fisher (Genentech)