

Fall 10-20-2015

Flexural fatigue behavior of rocking bioreactor films

Donald Buckley
GE Global Research, buckley@research.ge.com

K. Amanda Misner
GE Global Research

Terry Saunders
GE Global Research

Andrew Burns
GE Global Research

Hongyi Zhou
GE Global Research

See next page for additional authors

Follow this and additional works at: <http://dc.engconfintl.org/biopoly>

 Part of the [Materials Science and Engineering Commons](#)

Recommended Citation

Donald Buckley, K. Amanda Misner, Terry Saunders, Andrew Burns, Hongyi Zhou, Eugene Boden, Susan Burke, Ross Acucena, and Mark Roberts, "Flexural fatigue behavior of rocking bioreactor films" in "Single-Use Technologies: Bridging Polymer Science to Biotechnology Applications", Ekta Mahajan, Genentech, Inc., USA Gary Lye, University College London, UK Eds, ECI Symposium Series, (2015). <http://dc.engconfintl.org/biopoly/17>

This Conference Proceeding is brought to you for free and open access by the Proceedings at ECI Digital Archives. It has been accepted for inclusion in Single-Use Technologies: Bridging Polymer Science to Biotechnology Applications by an authorized administrator of ECI Digital Archives. For more information, please contact franco@bepress.com.

Authors

Donald Buckley, K. Amanda Misner, Terry Saunders, Andrew Burns, Hongyi Zhou, Eugene Boden, Susan Burke, Ross Acucena, and Mark Roberts



Flex Fatigue in Bioprocess Films

GE Global Research

Don Buckley, Amanda Misner,
Hongyi Zhou, Terry Saunders,
Andrew Burns, Gene Boden

GE Healthcare BioProcess

Mark Roberts, Sara Ullsten,
Susan Burke, Ross Acucena

October 20, 2015

Imagination at work.

GEHC Single-Use Disposable BioProcess Products



Xcellerex™ XDR, Wave™ & Xuri™ Bioreactors
(PL-01077/026, Bioclear 10 & 11)



ReadyToProcess™/ReadyCircuit™
(ReadyKleer)



HyClone & PAA Media/Reagents
(CX-514, ReadyKleer, Flexboy)

Multilayer Bioprocess Films:

Range of use conditions pose unique challenges across product lines

- Xcellerex: Pressure
- HyClone/PAA/RTP: Abrasion
- Wave: Flex Fatigue



WAVE Bioreactor™ Products



Rockers, Pumps & Controls



Sensor-Enabled
Disposable CellBags

Rocking action enables gas exchange, fluid mixing & suspension for batch, fed-batch & perfusion culture of a wide range of cell types.



Flex Fatigue Resistance

Motivation

- Flexural fatigue of folded film during rocking → Film stresses can progress to loss of bag integrity
- Standard testing (Gelbo, tensile, etc.) do not reproduce fatigue mechanism
- Rare, but highly impactful failure mode

Test Method

- High Rock Rate/Angle (25rpm/9°)
- One 50L bag per rocker for >30 days
- Does not necessarily test bags to failure
- Low throughput, variable test




Goal: Develop a more rapid & reproducible flex fatigue test for Wave™ Bioreactor Films



Understanding Bag Mechanics on Wave Rocker




50L Cellbag™, 24L water, 25 rpm, $\pm 10^\circ$



BioClear 10 50L WAVE CellBag

WAVE Rocker Testing

Analysis of real-world fold geometries on 50L WAVE CellBags



- Quantified static & rocking dimple geometries to understand film stresses
- Mechanism: Dimples roll, pivot & elongate during rocking generating local stresses in the film


Combined tensile + bending action on film results in fatigue



Flex Fatigue Test Method Evolution

Test Evolution

- Conventional flex tests are uncorrelated to rocker results
 - ASTM 2D, Gelbo, etc.
- Need greater constraint of ductile films to induce failure
- 3D “minibag” devised to reflect flex fatigue in Cellbags



GRC Flex Fatigue Test Approach:

- Reproducible dimple formation → fixed geometry
- More severe fold geometry → increased stress
- Rapid cycling of the fold (1Hz) → faster failure



Flex Test Progression:

Stages of Failure Seen in Wave Rocker Test and Flex Test

Rocker Test
(BioClear 10 film)



Flex Testing
(BioClear 10 film)



whitening

delamination
and/or crack


film break and leak

Flex Test Reproduces Sequence of Events in path to failure



Film Geometry


Neutral Plane: Film flexure causes tensile & compressive stresses, which are minimized near the neutral plane (y_0)



$$y_0 = \frac{1}{2} \times \frac{\sum_{i=1}^n E_i \times (h_i^2 - h_{i-1}^2)}{\sum_{i=1}^n E_i \times (h_i - h_{i-1})}$$




h - layer thickness
E - Young's modulus
 y_0 - neutral plane

BioClear 10/11
(300 μ m)





Flex Fatigue Failure Analysis – BioClear 10

External View



Cross-Section



Test Parameters

Mounting



- Positioning – fixed grip distance
- Preload (2 lb. nominal)
- Gripping (mechanical or hydraulic)

Indent

- Depth & radius of curvature mimic, but exceed rocker geometry

Flexing

- Frequency: 1 Hz → 24 hr. test for current prototype films
- Flexing waveform: sinusoidal (analogous to Wave)
- Flexing amplitude: Exceeds dimple flex amplitude on rocker




Set parameters to imitate rocker dimple,
and to generate reproducible data



Parallel Testing: multi-sample “4X” Tester


- 4 samples mounted/ Indented/ tested simultaneously
- Periodic dye test to detect through-holes
- Acquire fatigue data quickly to develop high level of statistical confidence



Research tool to assess flex fatigue resistance



4X Flex Fatigue Testing – Prototype Variants



Discriminates flex behavior of films to > 95% confidence



Summary

- A single film for use across a wide range of applications requires a balance of properties → need rapid assessment of prototype film performance
- WAVE is an extreme example of film flex fatigue (bulk liquid shipment, intermodal RTP shipment, STRs) → Films designed to exceed WAVE needs will exceed needs for other applications
- Developed a simple constrained flex test enabling discrimination of flex fatigue resistance in candidate films with high statistical confidence
- Cross-sectional microscopy elucidates failure sequences of different films, to guide decisions on structure optimization for flex fatigue
- Rapid screening of film candidates enables thorough rocker testing of promising candidates → reduce design cycle time & increase confidence



