

Fall 11-2-2015

Novel compact cell settlers for continuous perfusion bioreactor cultures of microbial (and mammalian) cells

Dhinakar Kompala

Sudhin Biopharma Company, kompala@sudhinbio.com

Brian Batt

Sudhin Biopharma Company

Follow this and additional works at: http://dc.engconfintl.org/biomanufact_ii



Part of the [Biomedical Engineering and Bioengineering Commons](#)

Recommended Citation

Batt, B.C., R.H. Davis, and D.S. Kompala, *Biotechnology Progress*, 6: 458-464, 1990 Searles, J.A., P.W. Todd, and D.S. Kompala, *Biotechnology Progress* 10: 198-206, 1994. Probstein, R.F. and D. Yung, US Patent 4151084, April 1979.

This Conference Proceeding is brought to you for free and open access by the Proceedings at ECI Digital Archives. It has been accepted for inclusion in Integrated Continuous Biomanufacturing II by an authorized administrator of ECI Digital Archives. For more information, please contact franco@bepress.com.

NOVEL COMPACT CELL SETTLERS FOR PERFUSION BIOREACTOR CULTURES OF MICROBIAL (AND MAMMALIAN) CELLS

Prem Singh S.D. Samuel, Cassandra A. Freeman,
Brian C. Batt, & Dhinakar S. Kompala

²Sudhin Biopharma Company, Superior, CO, USA

Kompala@sudhinbio.com

Summary of Abstract

As **secretory expression systems** have become well developed for some larger **microbial cells**, such as yeast *Saccharomyces cerevisiae* and *Pichia pastoris* cells, it is now possible to develop **high cell density continuous perfusion cultures** of *microbial yeast cells*.

Keywords: Compact cell settlers, Cell retention device, perfusion bioreactor, microbial yeast cells, selective removal of dead cells & cell debris, viable cell recycle.

Inclined Settler – a superior cell retention device

- Simple, passive, robust & powerful technology
- Uses no membrane barriers that can clog!
- Uses no moving parts that can break down!
- Proven at large scale for mammalian cell cultures
- Selective removal of dead cells and cell debris
- Complete recycle of live cells to bioreactor
- High *viable* cell density and productivity
- Perfusion bioreactors can be operated indefinitely
- Not demonstrated so far for microbial cells


Perfusion Bioreactor Culture of Yeast *P. pastoris* Cells using Inclined Settlers




Materials & Methods

- Cells: *Pichia pastoris* (x-33) size: 3 – 6 microns
- Bioreactor size: 1.5 liters, computer controlled
- Inclined Settler: CS-10 (from Dr. P. Brown, Biotechnology Solutions, Inc.) lamellar settler designed for 2 - 4 liter mammalian bioreactor
- Suggested max perfusion rate of ~8 liters/day for mammalian cell culture, i.e. 333 ml/hr
- Media: Glycerol, plus basal salts medium
- Temp: 28 °C pH: 5.6 Air flow: ~1 vvm


1st Perfusion Bioreactor - Pichia cells – 10% glycerol



Fraction of cells in settler outlet vs perfusion rate



Centrifuged samples from Settler outlet and Bioreactor



Lessons Learned from Lamellar Settler for the smaller yeast *Pichia* cells


- Inclined settler can separate dead cells and cell debris into the settler top outlet, while the larger live yeast cells are recycled back to the bioreactor.
- Very high cell densities (> 500 OD) are achieved and maintained for several weeks (interrupted by power outages in India shutting off pumps in weekends).
- CS-10 designed for up to 10 liters/day perfusion rate in mammalian cell cultures can be operated up to only 400 ml/day, before losing a significant fraction (>25%) of the smaller yeast cells.

Novel Compact Cell Settlers


- We have now built and patented some more efficiently scaled-up compact cell settlers.
- These novel cell settlers are smaller, need less footprint, and can scale up in true 3D.
- These novel cell settlers can be used for **microbial** and mammalian cell cultures.




Perfusion Bioreactor of *P. pastoris* cells with Compact Cell Settler as retention device




First Perfusion Bioreactor - Pichia cells using Compact Cell Settler as retention device



Size Distribution of Cells from Bioreactor and Settler Exit



Bioreactor & Settler Effluent Samples Centrifuged & Control



Conclusions

- Inclined settler technology can be successfully used to remove smaller dead cells / cell debris and achieve high cell densities for microbial perfusion cultures (as proven before for mammalian cultures).
- Rectilinear scale up of inclined settlers as lamellar settlers is inefficient for microbial cell cultures.
- More efficiently scaled up “compact cell settlers” can be used with microbial yeast cells, which are maintained at very high cell density (> 800 OD) for > 3 weeks so far in continuous perfusion bioreactors.

Key References

- Batt, B.C., R.H. Davis, and D.S. Kompala, *Biotechnology Progress*, 6: 458-464, 1990
- Searles, J.A., P.W. Todd, and D.S. Kompala, *Biotechnology Progress* 10: 198-206, 1994.

Acknowledgements

- National Science Foundation, Presidential Young Investigator award to DSK, \$500,000; 1988-1993
- National Science Foundation, Small Business Innovation Research Phase I grant \$150,000; 2015