

Fall 10-5-2015

A universal characterization method on viscous materials using depth sensing indentation

Abdul Shah

Jost Institute for Tribotechnology, akshah1@uclan.ac.uk

Nathalie Renevier

Jost Institute for Tribotechnology

Ian Sherrington

Jost Institute for Tribotechnology

Follow this and additional works at: http://dc.engconfintl.org/nanomechtest_v



Part of the [Materials Science and Engineering Commons](#)

Recommended Citation

1) BS EN ISO 14577-1:2015 2) Oliver, W.C., Pharr, G.M.:J. Mater. Res., 1992, 7, (06), pp. 1564–1583. 3) G. Feng, A.H.W.N.:J. Mater. Res., 2002, 17, (03), pp. 660 – 668. 4) Shah A, Renevier N, Sherrington I., UK Patent Application Number 1513480

This Abstract is brought to you for free and open access by the Proceedings at ECI Digital Archives. It has been accepted for inclusion in Nanomechanical Testing in Materials Research and Development V by an authorized administrator of ECI Digital Archives. For more information, please contact franco@bepress.com.

A UNIVERSAL CHARACTERIZATION METHOD ON VISCOUS MATERIALS USING DEPTH SENSING INDENTATION

Abdul Shah, Jost Institute for Tribotechnology, University of Central Lancashire, Preston, UK
akshah1@uclan.ac.uk

Nathalie Renevier, Jost Institute for Tribotechnology, University of Central Lancashire, Preston, UK
Ian Sherrington, Jost Institute for Tribotechnology, University of Central Lancashire, Preston, UK

Key Words: Nanoindentation; Method; Polymers; Creep; Plastic depth

Miniaturisation of devices and a paradigm shift toward using compliant material require small scale characterisation techniques such as nanoindentation [1]. Initial non-conformity of contact and delayed elasticity on the unloading curve are not currently taking into account in nanoindentation methods [2, 3], where the unloading curve is seen fully elastic. A different approach has been taken which considers actual localised deformation during nanoindentation, thus the proposed method [4] is able to acquire untainted elastic or viscoelastic response data. The method, validated for both viscous and non-viscous materials, takes into account the correction of the stiffness associated with the delayed elasticity at initial unloading and determines modulus with less variability even if testing conditions are non-quasi-static. The key step in the methodology is to account for initial conformity of the contact, the nose-out phenomena and delayed elasticity. Thus a full elastic point (FEP) is determined as seen in the figure 1. Traditionally to eliminate the initial delayed elastic response a multi-cycle approach was needed, whereas in our method any materials,




Figure 1 Load-displacement plot for PET

even viscous, can be tested under any test conditions in a single loading-unloading cycle. The algorithm is universally applicable. Experiments were conducted on six different viscous materials under single and multi cycle loading conditions to validate our method to existing ones [2–3]. Multi-cycle test on PET & PEN are reported in figure 2. It was found that except for rubber all the different materials studied by using just a single-cycle, our method determined values equally well as previous Oliver and Pharr method [2] with less variation. For multi-cycle tests our method is capable of producing results as good as the fourth-cycle of Feng's method [3], thus time and cost of experimentation can be reduced. So this method becomes appropriate as a standardised technique, and also for the characterisation of polymers which have been an issue in the past.




Figure 2 Modulus results against different methods using multi-cycle tests for a) PET & b) PEN

References

- 1) BS EN ISO 14577-1:2015
- 2) Oliver, W.C., Pharr, G.M.: J. Mater. Res., 1992, 7, (06), pp. 1564–1583.
- 3) G. Feng, A.H.W.N.: J. Mater. Res., 2002, 17, (03), pp. 660 – 668.
- 4) Shah A, Renevier N, Sherrington I., UK Patent Application Number 1513480.2.